

ENDO AND ENDO MICRO SURGERY KITS



ENDO AND ENDO MICRO SURGERY KITS

Инструменты для системы коффердам
Pinces porte-crampons a digue
Kofferdam Klammerzangen

Инструменты для системы коффердам
Pinzas de engrapar para dique de caucho
السد المطاط الملقط المشابك



Article# OI4-001
Anisworth



Article# OI4-002
Ivory



Article# OI4-003
Elastic
Seaparating & Ligature
Plier



Article# OI4-004
Elastic
Seaparating & Ligature
Plier



Article#
OI4-005
Brewer



Article#
OI4-006
Stockes



Article#
OI4-007
Palmer



Article#
OI4-008
Ivory
light weight



Article# OI4-009
Small
Rubber Dam Frame



Article# OI4-010
Medium
Rubber Dam Frame



Article# OI4-011
Large
Rubber Dam Frame

Инструменты для системы коффердам
Pincas porte-crampons a digue
Kofferdam Klammerzangen

Инструменти для системи коффердам
Pinzas de engrapar para dique de caucho
السد المطاط الملقط المشابيك



Article# OI4-012

On request the Content of the rubber dam kit can be changed with other instruments.

Content:

- 01 Anisworth
- 02 Brewer
- 03 Clamps Stand
- 04 Rubber Dam Frame L
- 05 Rubber Dam Frame M
- 06 Rubber Dam Box Lower Case
- 07 Rubber Dam Box Upper Case

- 08 Rubber Dam Clamp Fig. 212
- 09 Rubber Dam Clamp Fig. 211
- 10 Rubber Dam Clamp Fig. 8
- 11 Rubber Dam Clamp Fig. 8A
- 12 Rubber Dam Clamp Fig. 14
- 13 Rubber Dam Clamp Fig. 14A
- 14 Rubber Dam Clamp Fig. 0
- 15 Rubber Dam Clamp Fig. 00
- 16 Rubber Dam Clamp Fig. 2
- 17 Rubber Dam Clamp Fig. 2A
- 18 Rubber Dam Clamp Fig. 7



Article# OI4-013

On request the Content of the rubber dam kit can be changed with other instruments.

Content:

- 01 Ivory
- 02 Ivory light weight
- 03 Clamps Stand
- 04 Rubber Dam Frame M
- 05 Rubber Dam Frame L
- 06 Rubber Dam Box Lower Case
- 07 Rubber Dam Box Upper Case


- 08 Rubber Dam Clamp Fig. 209
- 09 Rubber Dam Clamp Fig. 212
- 10 Rubber Dam Clamp Fig. 211
- 11 Rubber Dam Clamp Fig. 8
- 12 Rubber Dam Clamp Fig. 8A
- 13 Rubber Dam Clamp Fig. 14
- 14 Rubber Dam Clamp Fig. 14A
- 15 Rubber Dam Clamp Fig. 0
- 16 Rubber Dam Clamp Fig. 00
- 17 Rubber Dam Clamp Fig. 2
- 18 Rubber Dam Clamp Fig. 2A
- 19 Rubber Dam Clamp Fig. 7



Article# OI4-014

Rubber Dam Set of 12 Pcs

- 01 Anisworth
- 02 Rubber Dam Frame
- 03 Ivory
- 04 Clamp
- 05 Clamp
- 06 Clamp
- 07 Clamp
- 08 Clamp
- 09 Clamp
- 10 Clamp
- 11 Clamp
- 12 Clamp



Article# OI4-015

Rubber Dum Template

Кламмера для коффердама
Crampons pour digue
Kofferdam Klammern

Кламмера для коффердама
Grapas para dique de caucho
المشابك السد المطاط



Article# OI4-016

Stainless Steel Gold Plated Rubber Dam Clamp Tray with 11 Different Clamps
Fig. 212, 211, 8, 8A, 14, 14A, 0, 00, 2, 2A, 7



Article# OI4-017

Stainless Steel Rubber Dam Clamp Tray with 12 Different Clamps
Fig. 209, 212, 211, 8, 8A, 14, 14A, 0, 00, 2, 2A, 7



Article# OI4-018

Anisworth 17cm, Ivory 17cm
Brewer 17.5cm, Rubber Dam Frame (Children)
Rubber Dam Frame (Adult), Rubber Dam Clamps 11 Pcs
Fig. 00, Fig.0, Fig.1, Fig.2A, Fig.3, Fig.7,
Fig.8, Fig.8A, Fig.9, Fig.212, Fig.14



Article# OI4-019

Ruber Dam Clamps 13 Pcs Black Color Tray
B4, 212L, 212R, W00, W2, W3, W7, W8A,
W8AD, 12A, 13A, 26N, 27 Long Black



Article# OI4-020

Stainless Steel Rubber Dam Clamp Tray with 8 Different Clamps
Fig. 1, 26, 29, 51, 60, 202, 205, 210



Article# OI4-021

Stainless Steel Rubber Dam Clamp Tray
(without clamps)

Rubber Dam Clamps Organizing Board



Article# OI4-022

Rubber Dam Clamp Organizing Board,
For 8 Clamps




Article# OI4-023

Rubber Dam Clamp Organizing Board,
For 12 Clamps

Кламмера для коффердама
 Crampons pour dique
 Kofferdam Klammern

Кламмера для коффердама
 Grapas para dique de caucho
 المشابك السد المطاط

All rubber dam clamps are made from surgical stainless steel.
 Matt finish of these Satin Steel Rubber Dam Clamps offset bright clinical lighting. Hand-crafted to absorb and dissipate light instead of reflecting it back. This section includes the Anterior Rubber Clamps.



Article# OI4-024
00
 Flat jaws and high bow for small lower canines and first premolars



Article# OI4-025
6
 Slightly offset jaws for labial caries on anterior teeth



Article# OI4-026
9
 Offset flat jaws and stiff bow for anterior teeth



Article# OI4-027
9S
 Offset flat serrated jaws for anterior teeth.




Article# OI4-028
212SA
 Slightly offset staggered jaws for odd shaped anterior.



Article# OI4-029
212
 Flat jaws for labial caries on anterior teeth.

Rubber Dam Clamps -Molar (General)

All rubber dam clamps are made from surgical stainless steel.
 Matt finish of these Satin Steel Rubber Dam Clamps offset bright clinical lighting. Hand-crafted to absorb and dissipate light instead of reflecting it back. This section includes the Molar (General) Rubber Clamps.



Article# OI4-030
8A
 Deeply festooned four-point jaws for molars



Article# OI4-031
W8A
 Slightly festooned jaws for molars.



Article# OI4-032
10
 Flat jaws designed for small upper left and lower right molars.



Article# OI4-033
11
 Flat jaws designed for small upper right & lower left molars



Article# OI4-034
14
 Deeply festooned jaws for partially erupted molars



Article# OI4-035
14A
 Deeply festooned jaws for partially erupted or irregularly shaped molars



Article# OI4-036
W14A
 Deeply festooned jaws for partially erupted or irregularly shaped molars



Article# OI4-037
18
 Flat jaws and stiff bow for molars



Article# OI4-038
24
 Flat jaws for upper right and lower left molars



Article# OI4-039
25
 Flat jaws for upper left and lower right molars



Article# OI4-040
27
 Flat jaws for large molars





Article# OI4-041
56
 Slightly festooned jaws for large molars



Article# OI4-042
56S
 Slightly festooned serrated jaws for large partially erupted molars

Кламмера для коффердама
Crampons pour digue
Kofferdam Klammern

Кламмера для коффердама
Grapas para dique de caucho
المشابك السد المطاط



Rubber Dam Clamps -Premolar

All rubber dam clamps are made from surgical stainless steel. Matt finish of these Satin Steel Rubber Dam Clamps offset bright clinical lighting. Hand-crafted to absorb and dissipate light instead of reflecting it back. This section includes the Premolar Rubber Clamps.



Article#
OI4-043
0

Flat jaws and high bow for longer premolars.



Article#
OI4-044
1

Slightly festooned jaws for upper premolars.



Article#
OI4-045
1A

Deeply festooned jaws for canines and premolars.



Article#
OI4-046
2

Flat jaws for premolars.



Article#
OI4-047
2A

Flat jaws for larger premolars.



Article#
OI4-048
27N

Slightly festooned jaws for premolars.



Article#
OI4-049
29

Flat jaws for lower premolars.



Article#
OI4-050
209

Small, flat jaws with stiff bow for lower premolars.



Article#
OI4-051
W2

Flat jaws for premolars.

Rubber Dam Clamps -Upper Molar

All rubber dam clamps are made from surgical stainless steel. Matt finish of these Satin Steel Rubber Dam Clamps offset bright clinical lighting. Hand-crafted to absorb and dissipate light instead of reflecting it back. This section includes the Upper Molar Rubber Clamps.



Article#
OI4-052
4

Flat narrow jaws for small upper molars.



Article#
OI4-053
W4

Flat narrow jaws for small upper molars.



Article#
OI4-054
5

Slightly festooned jaws for large upper molars.



Article#
OI4-055
8

Deeply festooned jaws for upper molars.



Article#
OI4-056
26N

Flat jaws for upper molars.



Article#
OI4-057
30

Offset jaws for upper right molars.



Article#
OI4-058
31

Offset jaws for upper left molars.



Article#
OI4-059
201

Slightly festooned jaws for upper molars.



Article#
OI4-060
205


Slightly festooned jaws for large upper molars.

Кламмера для коффердама
Crampons pour digue
Kofferdam Klammern

Кламмера для коффердама
Grapas para dique de caucho
المشابك السد المطاط

Rubber Dam Clamps -Lower Molar

All rubber dam clamps are made from surgical stainless steel. Matt finish of these Satin Steel Rubber Dam Clamps offset bright clinical lighting. Hand-crafted to absorb and dissipate light instead of reflecting it back. This section includes the Lower Molar Rubber Clamps.



Article#
OI4-061
3

Flat narrow jaws for small lower molars.



Article#
OI4-062
7

Flat jaws for lower molars



Article#
OI4-063
7A

Slightly reversed festooned jaws for large lower molars



Article#
OI4-064
12A

Serrated offset jaws for lower right molars



Article#
OI4-065
13A

Serrated offset jaws for lower left molars



Article#
OI4-066
26

Slightly reversed festooned jaws for small lower molars



Article#
OI4-067
28

Flat jaws for narrow molars



Article#
OI4-068
W3

Flat narrow jaws for small lower molars




Article#
OI4-069
W7

Flat jaws for lower molars.



Article# OI4-070
Plier For Clamp



Article# OI4-071
Enc Clamp Kit 11 Pcs

Кламмера для коффердама
Crampons pour digue
Kofferdam Klammern

Кламмера для коффердама
Grapas para dique de caucho
المشابك السد المطاط



Article# OI4-072 Fig. 15
Article# OI4-073 Fig. 18
Article# OI4-074 Fig. 24
Article# OI4-075 Fig. 25
Article# OI4-076 Fig. 26



Article# OI4-077 Fig. 28
Article# OI4-078 Fig. 30
Article# OI4-079 Fig. 31
Article# OI4-080 Fig. 51



Article# OI4-081 Fig. 138
Article# OI4-082 Fig. 139
Article# OI4-083 Fig. 200
Article# OI4-084 Fig. 201
Article# OI4-085 Fig. 202



Article# OI4-086 Fig. 203
Article# OI4-087 Fig. 204
Article# OI4-088 Fig. 205
Article# OI4-089 Fig. 3



Article# OI4-090 Fig. W3
Article# OI4-091 Fig. 4
Article# OI4-092 Fig. 56A
Article# OI4-093 Fig. 7
Article# OI4-094 Fig. 7A



Article# OI4-095 Fig. 8
Article# OI4-096 Fig. 8A
Article# OI4-097 Fig. W8A
Article# OI4-098 Fig. 10



Article# OI4-099 Fig. 11
Article# OI4-100 Fig. 12A
Article# OI4-101 Fig. 13A
Article# OI4-102 Fig. 14
Article# OI4-103 Fig. 14A



Article# OI4-104 Fig. 14AN
Article# OI4-105 Fig. 22
Article# OI4-106 Fig. 27
Article# OI4-107 Fig. 207



Article# OI4-108 Fig. 206
Article# OI4-109 Fig. 208
Article# OI4-110 Fig. 209
Article# OI4-111 Fig. 209
Article# OI4-112 Fig. 1



Article# OI4-113 Fig. 1A
Article# OI4-114 Fig. 2
Article# OI4-115 Fig. 2A
Article# OI4-116 Fig. W2A



Article# OI4-117 Fig. 0
Article# OI4-118 Fig. 00
Article# OI4-119 Fig. 210
Article# OI4-120 Fig. 211



Article# OI4-121 Fig. 212
Article# OI4-122 Fig. 6
Article# OI4-123 Fig. 9



Article# OI4-124 Fig. N06



Article# OI4-125 Fig. U67



Article# OI4-126 Fig. P1



Article# OI4-127 Fig. P2



Article# OI4-128 Fig. G-1



Article# OI4-129 Fig. G-2

Инструменты для удаления коронок
Arrache-couronnes
Teleskopkronenzangen

Инструменты для видалення коронок
Arranca-corona
مزيلات التاج



Article# OI4-130
Crown Spreaders



Article# OI4-131
Crown Spreaders
with Spring



Article# OI4-132
Wynman Crown Gripper
No Lock T.C 16 cm



Article# OI4-133
Crown Spreader



Article# OI4-134
Fig. 1



Article# OI4-135



Article# OI4-136
Crown Spreaders
Fig. 3



Article# OI4-137
Offset Crown Splitter
Remover Elevator



Article# OI4-138
Christenson crown remover
Straight, facial use anterior



Article# OI4-139
Christenson crown remover
Angled, for posterior use

Инструменты для удаления коронок
Arrache-couronnes
Teleskopkronenzangen

Инструменты для видалення коронок
Arranca-corona
مزيلات التاج



Article# OI4-140
16cm

Crown Plier with plastic inserts Straight
For ceramical and plastic crowns and bridges, with
plastic inserts (sterilization up to 180° C).



Article# OI4-141
14.5cm

Crown Plier
with plastic inserts
Angled



Article# OI4-142
14cm

Universal plier for the removal
of ceramic crowns and bridges



Article# OI4-143
14cm

Universal plier for the removal
of ceramic crowns and bridges



Article# OI4-144
Furrer
Copper Ring
Removing Plier.
15cm Curved



Article#
OI4-145 12.5 cm
OI4-146 15.0 cm
Bohm Pliers for telescope crowns.
Diamond points coated.



Article# OI4-147
Nerve-Canal Pliers
straight



Article# OI4-148
Nerve-Canal Pliers
Curved

Инструменты для удаления коронок
Arrache-couronnes
Teleskopkronenzangen

Инструменты для видалення коронок
Arranca-corona
مزيلات التاج



Article# OI4-149
Crown Remover
Set of 4 Pcs



Article# OI4-150
Crown Remover
Set of 4 Pcs



Article# OI4-151
Crown Remover
Set of 4 Pcs



Article# OI4-152
Crown Remover
Set of 4 Pcs



Article# OI4-153
Crown Remover
Set of 5 Pcs



Article# OI4-154
Crown Remover
Set of 6 Pcs



Article# OI4-155
Crown Spreader, Crown Remover
Special Instrument for crown spreading
after separating, simple and safe
operation.

Инструменты для удаления коронок
Arrache-couronnes
Teleskopkronenzangen

Инструменти для видалення коронок
Arranca-corona
مزيلات الناج



Fitting is possible on both ends

Fitting is possible on both ends

Article# OI4-156
Crown Remover,
For different level
of intensity



Single end Fitting

The inner spring can be remove, allow to control the stock manually

Article# OI4-157
Crown Remover,
For different level
of intensity



Fig. 1



Fig. 2



Fig. 3



Fig. 4



Fig. 5



Fig. 6



Article# OI4-158
Crown Remover
Set of 3 Pcs



A



B



Article# OI4-159
Set of 4 Pcs



Fig. 1



Fig. 2



Fig. 3

■ Матрицедержатели ■ Матрицеутримувачи
 ■ Porte Matrices ■ Porta Matrices
 ■ Matrizenspanner ■ مصفوفة الخدم



Article# OI4-160
Fig. 5
Ivory 8
5 mm

Article# OI4-161
Fig. 6
Ivory 8
6 mm

Article# OI4-162
Fig. 7
Ivory 8
7 mm

Article# OI4-163
Ivory 8 N

Article# OI4-164
Didner 8 P

Article# OI4-165
Didner 8 P
70 mm

Article# OI4-166
Matrix Ivory 8
6 mm

ADJUSTABLE



Article# OI4-167
Fig. 1
Tofflemire-Universal

Article# OI4-168
Fig. 2
Tofflemire-Junior

Article# OI4-169
Fig. 3
Tofflemire-Senior

Article# OI4-170
Nystrom
Left

Article# OI4-171
Nystrom
Right

Article# OI4-172
Narrow

Article# OI4-173
Wide


Article# OI4-174
Ivory
with shields

Article# OI4-175
Ivory
without shields



Article# OI4-176
Ivory

Article# OI4-177
Elliot




Article# OI4-178
Fig. 7

Article# OI4-179
Fig. 8

Article# OI4-180
Fig. 9

Article# OI4-181
Fig. 10

For Molars and Premolars



Article# OI4-182
Fig. 11

Article# OI4-183
Fig. 12

Article# OI4-184
Fig. 13
Tofflemire

Инструменты для имплантологии
Implantologie
Implantologie

Инструменти для імплантології
Implantologia
زراعة الأسنان



Article# OI4-185
0.8mm



Article# OI4-186
1mm



Article# OI4-187
1.2mm



Article# OI4-188
1.8mm



Article# OI4-189
Harvard Application
Forceps

Cement Teeth Filling Plugger Applicator
Carrier Root Canal Endodontic
Material 0.8/1.0/1.2/1.8/2.0/2.8/3.5mm



Article# OI4-190



Article# OI4-191



Article# OI4-192



Article# OI4-193



Article# OI4-194

■ Инструменты для амальгамы ■ Инструменти для амальгами
 ■ Instruments pour amalgame ■ Instrumentos para amalgama
 ■ Amalgamisinstrumente ■ الصكوك الملقم

Ø 2mm



Article# OI4-195
Fig. 10
16.5cm
English pattern

Ø 2mm



Article# OI4-196
Fig. 1
18.5cm
German pattern

Ø 2mm



Article# OI4-197
Fig. 1
15cm
Swiss pattern

Ø 2mm



Article# OI4-198
Fig. 1
15cm

Ø 2mm



Article# OI4-199
Fig. 2
15cm, 45°

Ø 2mm



Article# OI4-200
Fig. 3
14cm, 90°



Article# OI4-201
Chip Blower



Article# OI4-202



Article# OI4-203



Article# OI4-204



Article# OI4-205
Water Syringe

Инструменты для амальгамы Инструменти для амальгами
Instruments pour amalgame Instrumentos para amalgama
Amalgamisinstrumente الصكوك الملقم

AMALGAM CARRIERS

TIP SIZES FOR DISTAL OR SINGLE-END

1.5 mm

2.0 mm

TIP SIZES FOR DOUBLE-END

1.5 mm

2.8 mm

3.2 mm



Article#
OI4-206 - 1.5mm
OI4-207 - 2.0mm
OI4-208 - 2.0mm
OI4-209 - 3.0mm



Single-End



Article#
OI4-210
Mini/
Single-End
2.8 mm



Article#
OI4-212
Large
Single-End
2.0 mm



Article#
OI4-213
Regular
Distal



Article#
OI4-211
Regular
Single-End
1.5 mm



Article#
OI4-214
Mini
Distal



Article#
OI4-215
Mini/
Regular



Article#
OI4-216
Large/
Regular



Article#
OI4-217
Jumbo/
Regular



Article# OI4-230
Amalgam Well



Article# OI4-231
Stainless Steel
Amalgam Well
Medium



Article# OI4-232
Stainless Steel
Amalgam Well
large



Article# OI4-232A
Stainless Steel
Amalgam Well
large

TIP SIZES FOR DOUBLE-END



Article#
OI4-218 - 1.5mm-2.0mm
OI4-219 - 2.0mm-2.8mm
OI4-220 - 2.0mm-3.2mm



Article#
OI4-221 - 1.5mm-2.0mm
OI4-222 - 2.0mm-2.8mm
OI4-223 - 2.0mm-3.2mm




Article#
OI4-224 - 1.5mm-2.0mm
OI4-225 - 2.0mm-2.8mm
OI4-226 - 2.0mm-3.2mm



Article#
OI4-227 - 1.5mm-2.0mm
OI4-228 - 2.0mm-2.8mm
OI4-229 - 2.0mm-3.2mm


ENDOMICRO SURGERY KIT



©Association of Oral and Maxillofacial Surgeons of India



Root Canal Treatment



(a)



(b)



(c)



(d)



(e)



Article# OI4-246
Perry 9cm



Article# OI4-247
MTA Block




Article# OI4-248
Retractor




Article# OI4-251
Endomicro Surgery Kit with Cassette 19 Pcs

ENDOMICRO SURGERY SMALL KIT

Apicoectomy Treatment Stages:




Apicoectomy on upper front tooth



Micro apical surgery by SIE kit

Apical surgery is an option for the management of endodontically-treated tooth with persistent periapical lesions or symptom/sign. Several epidemiological studies have suggested that 33–60% of endodontically-treated teeth still presented the pictures of apical periodontitis. The possible causes may be persistent primary infection, secondary infection after endodontic therapy, vertical root fracture or cemental tears. The success rate of apical surgery was reported to range from 37% to 91%. Complete apical healing had been observed in 37%–96% of the endodontically-treated teeth after apical surgery.⁹ The wide range and inconsistency of the results may be attributed to the variation in treatment planning, surgical technique, methodology, and follow-up period.

Patient presented with moderate pain on percussion related to tooth 22. Upon radiographic and clinical examination we've found an apical lesion related to this tooth with normal PDL probing. The patient has done the RCT 15 years ago with the full arch prosthetic crowns that looks fine and healthy. So we discussed with the patient the different treatment plans and we suggest the micro apical surgery to be done for more conservative approach rather than removing the crown, post, Re-RCT, new post and new crown. The procedure was done in single visit compared to the multiple visits in the non surgical approach. The patient was very happy for her smile, so she wants to keep everything as it is, this could be accomplished by the surgical approach as well.



INSTRUMENTS SET FOR ENDODONTIC SURGERY

Article# OI4-263
Endomicro Surgery
Small Kit with
Cassette 14 Pcs

ENDOMICRO SURGERY SMALL KIT

INSTRUMENTS SET FOR ENDODONTIC SURGERY



(a)



(b)




(c)



(d)



(e)



Article# OI4-264 Explorer KX1, straight ø 0,45 mm



Article# OI4-265 Explorer KX2, angled ø 0,35 mm



Article# OI4-266 Micro plugger D1 ø 0,6 mm



Article# OI4-267 Micro plugger D2 ø 1,0 mm



Article# OI4-268 Filling spatula HEIDEMANN 1




Article# OI4-269 Excavator 33L



Article# OI4-270 Curette LUCAS 85



Article# OI4-271 Periosteal Elevator Prichard



Article# OI4-272 Mouth mirror with handle



Article# OI4-273 Micro mouth mirrors with handle (Ø 3 mm)



Article# OI4-278
Endomicro Surgery Small Kit
with Cassette 15 Pcs



Article# OI4-274
Needle holder Microvascular
with TC inserts





Article# OI4-275 Scalpel Handle Tip 16 cm Blade
Size 18, 20, 21, 22, 23, 24, 25, 60, 70



Article# OI4-276
De Backy Atraumatic 15cm Straight



Article# OI4-277
IRIS Scissors,
straight 9cm





(a)



(b)






(c)



(d)



(e)



Article# OI4-295 Micro Mirror

Article# OI4-296 Micro Mirror

Article# OI4-297 Micro Mirror

Article# OI4-298 Micro Mirror

MTAP

Article# OI4-299 MTA Carrier



Article# OI4-300

Potts-Smith Tissue Forceps 1x2 teeth 18 cm straight



Article# OI4-301

Potts-Smith Tissue Forceps 1x2 teeth 18cm curved



Article# OI4-302

Castroviejo Microsurgery Scissors-18cm



Article# OI4-303

Castroviejo Microsurgery Scissors-18cm



Article# OI4-304

Microsurgical Needle Holder Str-18cm



Article# OI4-305

Huey plier, diamond coated



Article# OI4-306

micro pliers with delicate flat tips, 90° and 45°



Article# OI4-307

Pin holding and nervebroaches plier

Article# OI4-308

Endo Micro Surgery
Complete Kit with Cassette 30 Pcs

